

AlcoMate® Premium



- ✓ **FDA 510(k) Cleared**
- ✓ **DOT / NHTSA Approved**
- ✓ **Meets US Coast Guard Requirements**
- ✓ **ISO 9002 Certified**

AlcoMate® Premium *Breath Alcohol Tester*

- Easy One Button Operation
- No Calibration Service Required
- Easy Replaceable PRISM Calibration Sensor
- Deep Lung Air Sampling
- Active Pressure Sensor
- Low Battery Indicator
- Internal Mouthpiece Storage Slot
- Self Diagnostic System
- Auto Power Off
- Range: 0 - 0.400% BAC

Innovative Drug Screening Solutions
Quality • Choices • Support

Learn More About
AlcoMate Employee Drug Testing Kits

 **EmployeeDrugTestingAce.com**
Trusted Proven Secure

www.DrugTestingAce.com


www.EDTAce.com

PREMIUM BREATHALYZER



Instant Technologies, Inc.



an  inverness medical company

AlcoMate® Premium

Premium

AlcoMate® Premium

Breath Alcohol Tester

The AlcoMate® Premium is our flagship breathalyzer that provides you with the most advanced, accurate and convenient commercial breath testing system available. **The Premium never requires calibration service.** With the unique Easy-Replace Sensor Panel, module replacement takes just a few seconds!



Target Users

Government Approvals
for Law Enforcement

Trusted for Use By

DOT (Department of Transportation
NHTSA (National Highway Traffic &
Safety Administration)

- Law Enforcement
- Medical Professionals
- School Administrators
- Corporate Employers
- Drug Rehabilitators
- Prosumers / General Consumers

The AlcoMate® Premium's PRISM™ Sensor Technology uses advanced next - generation components to improve upon prior oxide semiconductor technologies and introduces a level of accuracy, reliability and ease-of-use never before achieved among portable screening breathalyzers.

Note: The only legally valid measure of %BAC (percent Blood Alcohol Concentration) is obtained by blood testing or other evidential method of measurement. This device is for screening purposes only and should be used as a monitoring guide. This device must be used strictly in accordance with the instructions provided with the product.

Instructions

- 1 Press and hold power button.
The test count will be displayed.

- 2 Release power button and a countdown will begin. The display will show "BLOW" when ready.

- 3 Secure your lips around the mouthpiece and blow sharply and steadily into unit for about 4 seconds.



Power Button



Mouthpiece Slot

The results will be shown for about 15 seconds. The unit will then shut itself off.

Device Specifications:

Alcohol Sensor: Oxide Semiconductor Sensor
 Detection Range: 0.000% – 0.400% BAC
 Breath Sampling Time: 4 + seconds (Deep-Lung Air)
 Calibration: None (PRISM Technology)

Sensor Accuracy: + / - 0.01 at 0.10% BAC
 Display: 4 - Digit Decimal LCD (0.000%)
 Battery: Two 1.5V AA alkaline (included)
 Dimensions: 4.25 x 2.0 x 0.75 inches
 Casing: Silver Plastic, Light Grey Trim